
## ISCHEMIC STROKE SYSTEM



CEREBASE™DA guide sheath



LARGE BORE CATHETER



EMBOTRAP® || revascularization device



EMBOTRAP® II Revaso	EMBOTRAP® II Revascularization Device							
Product Code	Name	Recommended Vessel Diameter	Microcatheter Compatibility	Working Length	Device Length			
ET009521	5×21	1.5-5.0 mm	0.021"	21mm	194cm			
ET009533	5x33	1.5-5.0 mm	0.021"	33mm	195cm			

Large Bore Catheter						
Product Code	Distal ID	Distal OD	Proximal OD	Guide Compatibility	Usable Length	
IC71125UG	0.071"	0.081"	0.0825"	6F	125cm	
IC71132UG	0.071"	0.081"	0.0825"	6F	132cm	

CEREBASE™ DA Guide Sh	CEREBASE™ DA Guide Sheath						
Product Code	ID	OD (compatible short sheath)	Usable Length				
GS9070SD	0.090"	8F	70cm				
GS9080SD	0.090"	8F	80cm				
GS9090SD	0.090"	8F	90cm				
GS9095SD	0.090"	8F	95cm				



## To order contact customer service: Customer Service (800-255-2500 Option 1)

## Legal Manufacturer:

Medos International SARL Chemin-Blanc 38 CH-2400 Le Locle Switzerland

Neuravi Ltd. Block 3, Ballybrit Business Park Galway, H91 K5YD Ireland **Important information:** Prior to use, refer to the instructions for use supplied with this device for indications, contraindications, side effects, warnings and precautions.

Caution: US law restricts this device to sale by or on the order of a physician.