

“One size fits all.”

Great for tube socks,
not so great for
contact lens comfort.



extremeH2O[®] daily


Stop drop out rates and gain more referrals along the way. With three lens diameters for a perfect fit and long lasting comfort, you'll gain a contact wearer for life. Someone who loves their contacts, brings referrals and stays with the program.



Extreme H2O Daily is the ONLY one day disposable available in 13.6 14.2 14.8 . This enables you to successfully fit your small, medium and large cornea patients with the right fit, the first time - to improve comfort.

HVID Study Findings

Pacific University conducted a study on corneal diameter and findings indicate that 27% of contact lens fits fall outside the standard parameters and would benefit from a small or large diameter.



Material Benefits

Extreme H2O Daily is made with a non-ionic, ultra hydrating material that retains 99% of water content on the eye through the entire wearing time. Lens comfort, movement and fit are stable throughout the day.

EXTREME H2O DAILY PRODUCT SPECIFICATIONS

Material / Handling Tint	41% hioxifilcon A, 59% water / Light blue handling tint		
Polymer Type	GMA/HEMA copolymer, group 2, high water, non-ionic		
O ₂ Permeability	28 Dk (Fatt units @ 35°C)		
Packaging	Blister pack - 30 per box and/or 90 per box		
Wear Indication	Daily replacement suggested		
Center Thickness	0.070 mm @ -3.00		
Diameter	13.6	14.2	14.8
Base Curve	Median / 8.2	Median / 8.5	Median / 8.7
Powers	+0.50 to +6.00 -0.50 to -10.00	+0.50 to +6.00 -0.50 to -10.00	+0.50 to +6.00 -0.50 to -10.00
<i>.50D Steps after -6.00</i>			

Georgia: 2775 Premiere Parkway, Suite 600, Duluth, GA 30097 | 800.241.9312
 Florida: 7575 Commerce Court, Sarasota, FL 34243 | 877.336.2482
www.xcelspecialtycontacts.com

