

Size matters.



extremeH2O[®] weekly

Stop drop out rates and gain more referrals along the way. With three lens diameters for a perfect fit and long lasting comfort, you'll gain a contact wearer for life. Someone who loves their contacts, brings referrals, and stays with the program. It's simple, without Extreme H2O you're leaving money on the table.



Extreme H2O Weekly is the ONLY disposable available in three diameters: 13.6, 14.2 and 14.8. This enables you to successfully fit your small, medium and large cornea patients with the right fit, the first time, to ensure comfort.

13.6 14.2 14.8



X-CEL


SPECIALTY CONTACTS

Contact X-Cel Specialty Contacts to learn more.

Main: 877.336.2482 | www.xcelspecialtycontacts.com

HVID Study Findings

Pacific University conducted a study on corneal diameter and findings indicate that 27% of contact lens fits fall outside the standard parameters and would benefit from a small or large diameter.



Material Benefits

Extreme H2O Weekly is made with a non-ionic, ultra-hydrating material that retains 99% of water content on the eye through the entire wearing time. Lens comfort, movement and fit are stable throughout the day.

EXTREME H2O WEEKLY PRODUCT SPECIFICATIONS

Material / Handling Tint	41% hioxifilcon A, 59% water / light blue handling tint			
Polymer Type	GMA/HEMA copolymer, group 2, high water, non-ionic			
O ₂ Permeability	28 Dk (Fatt units @ 35°C)			
Packaging	Blister pack - 12 per box			
Wear Indication	Daily wear / weekly replacement suggested			
Center Thickness	0.110 mm @ -3.00			
Diameter	13.6	14.2		14.8
Base Curve	Median / 8.2	Median / 8.5	Steep / 8.3	Median / 8.7
Powers <i>.50D Steps after -6.00</i>	+0.50 to +6.00 -0.50 to -10.00	+0.50 to +6.00 -0.50 to -10.00	-0.50 to -10.00	+0.50 to +6.00 -0.50 to -10.00

Contact X-Cel Specialty Contacts to learn more.
Main: 877.336.2482 | www.xcelspecialtycontacts.com

