Simplified Post Surgical Fitting!





Advanced Design. Simple Success.

Intense research and development at Metro Optics has yielded a simple and effective approach to addressing the needs of the post surgical patient. Working with eye care professionals throughout the United States Metro Optics has created a fitting system for post surgery that can be easily used in any eye care practice.

RSa

Versatile enough to fit all post surgical conditions, this lens easily and consistently provides patients with excellent comfort and visual acuity while fitting within the corneal surface. A combination of spherical and aspheric curves create alignment with the peripheral cornea and slight vaulting of the surgically altered central cornea. Standard diameters are used for PKP patients, and larger diameters are used for post refractive patients.

Toric parameters can also be incorporated into this lens design to address patients with significant astigmatism.

Metro Optics... Always Ahead of the Curve

The Metro Optics Advantage

- Simple Warranty/Return Policy
- Free Standard Shipping
- Free Plasma Treatment on all GP Lenses
- Free Dot on GP Lenses
- Free Minus Lenticular

RSb

This lens features all the versatility of the RSa in consistently providing patients with excellent comfort and visual acuity. It also features secondary base curve geometry for improved edge clearance on post refractive patients having excessively steep peripheral corneas.

A larger than standard diameter is used to aid lens centering.

Toric parameters can also be incorporated into this lens design to address patients with significant astigmatism.

Signature Series

