

VISULENS 500 from ZEISS Automated Lensmeter

Make your day-to-day work easier



The easiest way to measure lenses



Make your day-to-day work easier and ensure precise measurement with an automated lensmeter from ZEISS

Using the VISULENS® 500 from ZEISS you can speed up the refraction and prescription processes.

The ZEISS VISULENS 500 is equipped with a green measurement LED. It allows you to precisely determine the power of the entire set of different lens materials and designs without the need to adjust the Abbe number. This benefit will add speed as well as reliability to your examination process. In addition, ZEISS VISULENS 500 features an advanced UV transmission measurement system using four different wavelengths (365, 375, 395, 405 nm). The result: increased confidence for your patient and more upselling support for you.

Make your life easier and increase patient satisfaction at the same time

The ZEISS VISULENS 500 is a straightforward, easy-to-use automated lensmeter designed to speed up your workflow. The instrument helps you to spend more time on client and patient consultation, or simply to increase your throughput. It will support you by checking correct lens centration.

Highlights at a glance



'Going green' allows easy and precise measurement of entire sets of different lens materials and indices without the need to adjust the Abbe number. Simply take a measurement and have confidence in the results. Analyze soft and hard contact lenses using the adapter accessories.



Comprehensively analyze lenses for their UV protective properties by measuring UV transmission simultaneously at four different wavelengths between 365-405 nm.



The ZEISS VISULENS 500 provides a 3-pin lens holder for the secure and reliable mounting of any kind of lens. A precise pen-type system ensures easy operation and clear marking of optical center and cylinder axis.

Connect your device to a local network, use a patient work list on a tablet and transfer measurements to your EMR/PMS or FORUM® from ZEISS with the optional VISUCONNECT® 500 from ZEISS.



Further benefits

- Holder for frames adjustable for a wide range of frames (straight, loop, rimless, nylon)
- Optional contact lens adapter with precise
 XY positioning
- Integrated thermo printer
- Tiltable TFT display for operator comfort
- PD measurement
- Serial data output via RS232 interface
- Progressive/bifocal lens mode



Technical Data

VISULENS 500 from ZEISS

6 different operating modes are available

- Standard mode for single vision, multifocal and prism lenses
- Tinted lens mode
- Progressive lens mode
- Soft contact lens mode
- Hard contact lens mode
- UV transmission mode

Measurement wavelength 546 nm (e line)

Measuring range

Sphere power	-25 D to +25 D
Cylinder power	0 D to ±10.00 D
Cylinder axis	0° to 180°
Addition value	0 D to +10 D
Prism value	0 Δ to 10 Δ
Pupil distance	0 mm to 81.8 mm

Measurement display

Diopters	0.01/0.125/0.25 D
Cylinder	Mix, +, -
Cylinder axis	1°
Prisms	0.01/0.125/0.25 Δ Cartesian/polar/ decentering
Display	Tiltable 320 x 240 color TFT-LCD
Data interface	RS232
Power supply	100 V to 240 V~, 50/60 Hz, 40 VA
Dimensions (W x D x H)	210 mm x 270 mm x 417 mm
Weight	6 kg

UV-Transmission measurement

Wavelength 365 nm, 375 nm, 395 nm, 405 nm	
---	--



ZEISS VISULENS 500 delivers accuracy and reliability right from day one

C 6 0297

Manufactured for:

Carl Zeiss Meditec AG Goeschwitzer Strasse 51-52 07745 Jena Germany



Carl Zeiss Meditec, Inc.

5160 Hacienda Drive Dublin, CA 94568 USA www.zeiss.com/us/med

from the current status of approval of the product or service offering in your country. Please contact our regional representatives for more information. Subject to change in design and scope of delivery and as a result of ongoing technical development. VISULENS and VISUCONNECT are either trademarks or registered trademarks of United States Edition: Only for sale in selected countries. The contents of the reference guide may differ