



VISUREF 100 from ZEISS
Autorefractor / Keratometer

Premium vision starts with premium diagnostics



We make it visible.

A valuable addition to your practice



What you face every day

When taking care of patients, every minute is valuable. There is no time for dealing with complex systems. Still, you need diagnostic data that is reliable and accurate, no matter what. Especially for the very first step in vision care: the refraction/keratometric evaluation.

How we can help

VISUREF® 100 from ZEISS – a straightforward, easy-to-use diagnostic instrument provides accurate and reliable data for your daily practice routines.

Easy to set up and use, the ZEISS VISUREF 100 can be operated after only a short training period. Intuitive handling features let you conduct the most important measurements with one device, including contact lens base curve, pupil and iris diameter evaluations. Refraction measurements of patients with implanted IOLs are also possible. Additionally, the automatic fogging incorporates an infinity scene to reduce the effects of accommodation.

As such, the ZEISS VISUREF 100 quickly becomes a valuable part of your practice environment, adding more precision and reliability to your vision care – right from the start.

Highlights at a glance

Motorized chinrest

Conveniently operated by two buttons next to the joystick, the motorized chinrest lets you quickly and comfortably position the patient.

Tiltable, 6.4 inch TFT display

Diagnostic information is clearly displayed in excellent contrast and resolution on the adjustable TFT monitor. Featuring a tilt range of 20°, it can be adjusted for individual user viewing preferences.

Quick release brake

Within easy reach of the joystick, the quick release brake can be used to disable the XY movement of the instrument base. It helps you to comfortably position the instrument and increase the ease of use when the patient is seated.

Connectivity with ZEISS VISUCONNECT 500

Connect your device to a local network, use a patient work list on a tablet and transfer measurements to your EMR/PMS or FORUM® from ZEISS with the optional VISUCONNECT® 500 from ZEISS.

Key benefits for you and your patients

- Accurate and reliable data for high quality patient care
- Convenience in measurement and handling
- Transparent information and clear visibility for data analysis
- Patient comfort through ergonomic design
- The optional ZEISS VISUCONNECT 500 provides network connectivity, direct DICOM and a well-defined EMR interface



Technical data

VISUREF 100 from ZEISS

4 different operation modes

- Refraction and keratometry combined (RK)
- Refraction (REF)
- Keratometry (KER)
- Contact lens base curve (CLBC)

Refraction

Sphere (SPH)	-25 D to + 22 D in 0.12/0.25 D increments
Cylinder (CYL)	0 D to ± 10 D in 0.12/0.25 D increments
Axis (AX)	0°–180° in 1° increments
Vertex distance (VD)	0.0, 10.0, 12.0, 13.5, 15.0 mm
Pupil distance (PD)	10–85 mm
Min. pupil diameter	2 mm

Keratometry

Corneal curvature	5.0–10.2 mm in 0.01 mm increments
Corneal refraction	33 D to 67.5 D in 0.12/0.25 D increments
Corneal astigmatism	0 D to -15 D in 0.12/0.25 D increments
Axis	0°–180° in 1° increments

Corneal diameter	2.0–12.0 mm in 0.1 mm increments
Chinrest movement	65 mm, motorized
Printer	Thermal, integrated
Connectivity	1 x RS232, 1 x USB, 1 x ext. video
Power supply	AC 100V–240 V, 50/60 Hz
Dimensions, weight	275 mm x 525 mm x 450 mm (WDH), 18 kg



ZEISS VISUREF 100 delivers accuracy and reliability right from the start



0297



Carl Zeiss Meditec AG
Goeschwitzer Strasse 51–52
07745 Jena
Germany
www.zeiss.com/visuref100
www.zeiss.com/med/contacts

EN_33_022_00031 Printed in Germany. CZ-11/2015 International Edition: Only for sale in selected countries.
The contents of the brochure may differ from the current status of approval of the product in your country. Please contact your regional representative for more information.
Subject to change in design and scope of delivery and as a result of ongoing technical development. VISUREF and VISUCONNECT are either trademarks or registered trademarks of Carl Zeiss Meditec AG or other companies of the Carl Zeiss Group in Germany and/or other countries.
© Carl Zeiss Meditec AG, 2015. All rights reserved.