## Introduction

The hippocampus (Hc), a fundamental neural substrate of memory, is composed of cytoarchitectonically distinct subfields: dentate gyrus (DG), cornu ammonis sectors (CA1-3) and subiculum [1], thought to support distinct memory aspects [4].

Evidence of differential maturation across the Hc subfields [2] suggests that Hc subfields maturation underlies memory development [3,4]. This possibility can be tested with longitudinal structural MRI studies, with reliable measurement of Hc subfield volumes, such as that achieved by manual demarcation.

To interpret differential development of these regions, error variance must be the same across subfields. Yet, Hc subfields structure can introduce potentially differential source of errors due to their distinct morphological and volumetric properties. Importantly, the errors may show age related bias in developmental samples due to the factors such as differences in movement or head size across visits. Manual demarcation additionally may exhibit a degree of variability over time and between raters, leading to inconsistency of Hc subfields measurement between time points[5].

Here, we assessed the test-retest consistency of Hc subfields volume measures obtained with manual demarcation of ultra-high Hc MR images acquired at different time points. We established a high longitudinal consistency of Hc subfield volumes using our methods, a prerequisite for interpreting meaningful longitudinal changes.

## **Participants**

Sample 1: two MRIs, one-month delay (n=28, ages 7-20 years, mean=12.64, SD=3.35)

**Sample 2:** two MRIs, two-year delay (n=28, ages 8-17 years, mean=11.90, SD=2.78)

### Brain Imaging

High-resolution structural MR images (T2 PD TSE, 0.4x0.4x2.0 mm<sup>3</sup>) were collected with a 3T Scanner at WSU

## Hippocampal Subfields Volume

Manual demarcation of the subfields were conducted by raters with high inter-rater reliability (ICC(2) > 0.85, ICC(3) > 0.90)

## Longitudinal Protocol


A randomized procedure was implemented to blind raters to the information about ID and timepoints. To make consistent decision across two timepoints, the brain pairs were traced simultaneously and slice by slice

## **Statistical Analyses**

Test-retest Reliability Using Intraclass Correlation Coefficient: Two way mixed-effect single ANOVA model that assumes non-independent observation (ICC3,1; Shrout and Fleiss, 1979 [7]) was conducted for Hc subfield volumes across two visits

**Bias Evaluation**: To determine if measurement error correlates with size of Hc subfields we implemented Bland-Altman Plots. Pearson correlation and one sample t-test were conducted to assess the measurement bias. Bar-graphs and boxplots were plotted to visualize the ICC(3) measures and sensitivity of the method in capturing individual differences, respectively. Error bars represent 95% confidence intervals.

## Methods



# **Test-Retest Consistency of Hippocampal Subfield Volume Measures: Implications for Longitudinal Developmental Studies**

Roya Homayouni, Qijing Yu, Sruthi Ramesh, Ana Daugherty, Noa Ofen Institute of Gerontology & Department of Psychology, Wayne State University

> These findings support excellent consistency between Hc subfield volume measures assessed by our manual demarcation protocol and MRI acquisition parameters.

**Over two years**, we found excellent consistency between Hc subfields volume measures, yet, compared to over one month, there was larger individual variability in DG-CA3 volumes, indicating high sensitivity of our method



Adolescents (age 13-20 years)

Subiculum: 2027.07 [818.93-3257.13]

- representations: CA2 joins the circuit. Trends Neurosci 34:526–535. doi:
- Lavenex, P., Lavenex P. B., (2013). Building hippocampal circuits to learn and remember: Insights into the development of human memory. Behav
- 3. Daugherty, A. M., Flinn, R. W., Ofen, N. (2017). Age-related differences in CA3-dentate gyrus volume uniquely linked to improvement in associative memory from childhood to adulthood. NeuroImage. 153: 75-85.
- Bergner, M., Newcombe. N.S. (2018). Hippocampal maturation drives memory from generalization to specificity. Trends Cogn. Sci., 22 (2018), рр. 676-686
- 5. Yushkevich P.A. Pluta J.B. Wang H. Xie L. Ding S.L. Gertje E.C. Mancuso L. Kliot D. Das S.R. Wolk D.A. (2015). Automated volumetry and regional thickness analysis of hippocampal subfields and medial temporal cortical structures in mild cognitive impairment. Hum. Brain Mapp; 36: 258-287
- 6. Shrout PE, Fleiss JL (1979): Intraclass correlations: Uses in assessing rater reliability. Psychol Bull 86:420–428